


Mass Spectrum



Sample Information

User : Gary
 Month-Day Processed : 07/20/16
 Time Processed : 10:32:50
 Injection Volume : 0.2
 Sample Name : C-Peptide (3-33), human
 Cat No. : RPI1116
 Lot No. : PE7711607
 Theoretical MW : 3020.26
 Observed MW : 3019.90

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250°C
 Block Temp : 200°C

Interface Bias : +4.5kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2ml/min
 B.conc : 50%H2O/50%MeOH