

CERTIFICATE OF ANALYSIS

Product Information

Product Name: FGF-basic (146 aa), Human

Cat. No.: Z03166-10; Z03166-50; Z03166-1

Lot No.: P50121412

Size: 10 µg/vial; 50 µg/vial; 1 mg/vial

Storage Condition: Lyophilized recombinant human Fibroblast Growth Factor-basic (146 aa) (rhFGF-basic) remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, rhFGF-basic remains stable up to 2 weeks at 4°C or up to 3 months at -20°C. For long term storage it is recommended that a carrier protein (example 0.1% BSA) can be added. Avoid repeated freeze-thaw cycles.

Quality Control

Test Items	Specifications	Results
Physical Appearance	White Freeze-dried Powder	Conformed
Purity (SDS-PAGE)	> 95%	Conformed
Biological Activity (ED ₅₀)	< 0.25 ng/mL	Conformed
Molecular Weight	16.4 kDa	Conformed
Bacterial Endotoxin	< 0.2 EU/µg	Conformed

Caution:

For laboratory use only. Not for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by: Joe Wang Date: 12/22/2014

Protein Products Department