


Mass Spectrum



Sample Information		Interface :ESI	Interface Bias :+4.5kV
User :	Gary	Nebulizing Gas Flow :1.5L/min	Drying Gas Flow :5 L/min
Month-Day Processed :	05/11/17	CDL Temp :250°C	T.Flow :0.2ml/min
Time Processed :	10:17:47	Block Temp :200°C	B.conc :50%H2O/50%MeOH
Injection Volume :	0.2		
Sample Name :	Glucagon-Like Peptide (GLP) I (7-37)		
Cat No. :	RP12738		
Lot No. :	PE9371705		
Theoretical MW :	3355.68		
Observed MW :	3356.40		