


Mass Spectrum



Sample Information

User : Gary
 Month-Day Processed : 06/30/17
 Time Processed : 14:32:17
 Injection Volume : 0.5
 Sample Name : Adrenomedullin (AM) (1-52), human
 Cat No. : RP11288
 Lot No. : PE2341706
 Theoretical MW : 6028.90
 Observed MW : 6029.10

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250 °C
 Block Temp : 200 °C

Interface Bias : +4.5kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2ml/min
 B.conc : 50%H2O/50%MeOH