


# Mass Spectrum



**Sample Information**  
Month-Day Processed : 02/03/21  
Time Processed : 2:17:17 PM  
Injection Volume : 0.3  
Sample Name : LCMV gp33-41  
Cat No. : RP20257  
Lot No. : E3755421B  
Theoretical MW : 1044.25  
Observed MW : 1044.0

**Interface** : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

**Equipment** : GK11010007  
Interface Bias : +3.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50%H2O/50%MeOH