


Mass Spectrum



Sample Information

Month-Day Processed : 11/07/20
Time Processed : 03:12:12 PM
Injection Volume : 0.3
Sample Name : C-Type Natriuretic Peptide (CNP) (1-22), human
Cat No. : RP111110
Lot No. : E8654820K
Theoretical MW : 2197.60
Observed MW : 2196.9

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH