


Mass Spectrum



Sample Information
Month-Day Processed : 05/24/21
Time Processed : 4:59:24 PM
Injection Volume : 0.4
Sample Name : mTRP-2(180-188)
Cat No. : RP20250
Lot No. : E6784021E
Theoretical MW : 1175.34
Observed MW : 1174.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH