


Mass Spectrum



Sample Information
Month-Day Processed : 10/31/20
Time Processed : 21:05:44
Injection Volume : 0.4
Sample Name : Exendin-4
Cat No. : RP10874
Lot No. : E3668520J
Theoretical MW : 4186.66
Observed MW : 4185.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH