


Mass Spectrum



Sample Information

User : Gary
 Month-Day Processed : 01/09/20
 Time Processed : 19:09:35
 Injection Volume : 0.3
 Sample Name : LL - 37, scrambled
 Cat No. : RP20332
 Lot No. : E5295320A
 Theoretical MW : 4493.30
 Observed MW : 4493.0

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment : GK11010007
 Interface Bias : +4.5kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2ml/min
 B.conc : 50%H2O/50%MeOH