


Mass Spectrum



Sample Information
 Month-Day Processed : 08/20/21
 Time Processed : 7:49:05
 Injection Volume : 0.2
 Sample Name : Neuropeptide Y, human
 Cat. No. : RP10492
 Lot. No. : E6321821G
 Theoretical MW : 4271.68
 Observed MW : 4272.0

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment : GK11010007
 Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50% H₂O/50% MeOH