


Mass Spectrum



Sample Information
Month-Day Processed : 11/24/20
Time Processed : 03:35:32 PM
Injection Volume : 0.4
Sample Name : MDYKDHGDYKDHIDYKDDDDK peptide
Cat No : RP21087
Lot No : E4648820K
Theoretical MW : 2861.87
Observed MW : 2861.5

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH