


Mass Spectrum



Sample Information

Month-Day Processed : 12/22/20
Time Processed : 10:28:22 PM
Injection Volume : 0.4
Sample Name : c-Myc Peptide
Cat No : RP11731
Sample ID : E4919820L
Theoretical MW : 1203.30
Observed MW : 1203.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH