


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 04/08/20  
Time Processed : 07:48:46 AM  
Injection Volume : 0.3  
Sample Name : Hexa-His  
Cat No. : RP11737  
Lot No. : E0401920D  
Theoretical MW : 840.85  
Observed MW : 840.6

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200  
Equipment : GK11010007  
Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH