


Mass Spectrum



Sample Information
Month-Day/Year Processed : 04/24/21
Time Processed : 11:24:16 AM
Injection Volume : 0.4
Sample Name : TAT peptide
Cat No. : RP20256
Lot No. : E7387821D
Theoretical MW : 1621.91
Observed MW : 1621.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH