


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 01/23/18
Time Processed : 19:47:33
Injection Volume : 0.2
Sample Name : RGD
Cat No. : RP20297
Lot No. : PE1631801
Theoretical MW : 1025.06
Observed MW : 1025.70

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH