

Sample Name :Flagelin22  
 Cat No. :RP19986  
 Lot No. :E2668221E  
 Time Processed :23:40:38  
 Month-Day-Year Processed :05/10/2021

Pump A : 0.065% trifluoroacetic in 100% water (v/v)  
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)  
 Total Flow:1 ml/min  
 Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm  
 Equipment: ZJ19010324

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	11.724	7164	1007	0.168
2	11.967	4059807	599816	95.051
3	12.285	155411	12623	3.639
4	12.694	28443	3091	0.666
5	13.109	20369	1760	0.477
Total		4271194	618296	100.000