


Mass Spectrum



Sample Information

Month-Day Processed : 10/11/20
Time Processed : 08:56:22 PM
Injection Volume : 0.3
Sample Name : PAR-4 Agonist Peptide, amide
Cat No. : RP11529
Lot No. : E2277320J
Theoretical MW : 680.80
Observed MW : 680.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH