


# Mass Spectrum



**Sample Information**  
 Acquired by : Gary  
 Month-Day Processed : 02/24/18  
 Time Processed : 16:43:08  
 Injection Volume : 0.4  
 Sample Name : Exendin-4  
 Cat No.: RP10874  
 Lot No.: PE4031802  
 Theoretical MW : 4186.66  
 Observed MW : 4186.80

**Interface** : ESI  
 Nebulizing Gas Flow : 1.5L/min  
 CDL Temp : 250°C  
 Block Temp : 200°C

**Interface Bias**  
 Drying Gas Flow : +4.5 kV  
 T.Flow : 5 L/min  
 B.conc : 0.2 ml/min

: 50% H<sub>2</sub>O/50% MeOH