


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 07/10/20  
Time Processed : 12:22:11 PM  
Injection Volume : 0.3  
Sample Name : MDYKDHGDYKDHIDYKDDDDDK peptide  
Cat No. : RP21087  
Lot No. : E6127720G  
Theoretical MW : 2861.87  
Observed MW : 2861.6

Equipment : GK11010007  
Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200