


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 08/29/19
Time Processed : 07:18:35 PM
Injection Volume : 0.3
Sample Name : Hexa-His
Cat. No. : RP11737
Lot. No. : E0016819H
Theoretical MW : 840.85
Observed MW : 840.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H_2O /50% MeOH