


# Mass Spectrum



**Sample Information**  
 Month-Day Processed : 05/31/21  
 Time Processed : 9:26:40  
 Injection Volume : 0.2  
 Sample Name : Neuropeptide Y, human  
 Cat. No. : RP10492  
 Lot. No. : E1820221E  
 Theoretical MW : 4271.68  
 Observed MW : 4272.6

**Interface** : ESI  
 Nebulizing Gas Flow : 1.5L/min  
 CDL Temp : 250  
 Block Temp : 200

**Equipment** : GK11010007  
 Interface Bias : +4.5 kV  
 Drying Gas Flow : 5 L/min  
 T.Flow : 0.2 ml/min  
 B.conc : 50% H<sub>2</sub>O/50% MeOH