


Mass Spectrum



Sample Information

Acquired by: Gary
 Month-Day Processed: 11/16/19
 Time Processed: 08:39:09 AM
 Injection Volume: 0.3
 Sample Name: Parathyroid Hormone (PTH) (1-34), Human
 Cat. No.: RP01001
 Lot. No.: E5737619K
 Theoretical MW: 4117.74
 Observed MW: 4117.0

Interface: ESI
 Nebulizing Gas Flow: 1.5L/min
 CDL Temp: 250
 Block Temp: 200

Equipment: GK11010007
 Interface Bias: +4.5 kV
 Drying Gas Flow: 5 L/min
 T.Flow: 0.2 ml/min
 B.conc: 50% H2O/50% MeOH