


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 03/06/17  
Time Processed : 15:19:15  
Injection Volume : 0.2  
Sample Name : TAT peptide  
Cat No. : RP20256  
Lot No. : PE7271702  
Theoretical MW : 1621.91  
Observed MW : 1622.20

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH