


# Mass Spectrum



## Sample Information

|                       |   |                       |            |
|-----------------------|---|-----------------------|------------|
| Month-Day Processed : | 04/30/21                                    | Equipment :           | GK11010007 |
| Time Processed :      | 7:21:36                                     | Interface :           | ESI        |
| Injection Volume :    | 0.4   | Nebulizing Gas Flow : | 1.5L/min   |
| Sample Name :         | Calcitonin Gene Related Peptide (CGRP), rat | CDL Temp :            | 250        |
| Cat. No. :            | RP11095                                     | Block Temp :          | 200        |
| Lot. No. :            | E9406121D                                   |                       |            |
| Theoretical MW :      | 3806.25                                     |                       |            |
| Observed MW :         | 3806.0                                      |                       |            |