


# Mass Spectrum



Sample Information  
Month-Day Processed : 03/14/21  
Time Processed : 12:18:08 AM  
Injection Volume : 0.5  
Sample Name : CEF20  
Cat. No. : RP11781  
Lot. No. : E8215221C  
Theoretical MW : 943.18  
Observed MW : 942.4

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200  
Equipment : GK11010007  
Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50%H2O/50%MeOH