


Mass Spectrum



Sample Information
 Month-Day Processed : 07/02/21
 Time Processed : 11:51:59
 Injection Volume : 0.2
 Sample Name : Calcitonin Gene Related Peptide (CGRP), rat
 Cat. No. : RP11095
 Lot. No. : E8359521F
 Theoretical MW : 3806.25
 Observed MW : 3806.0

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment : GK11010007
 Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50% H2O/50% MeOH