

Sample Name :Gastric Inhibitory Peptide (GIP), human  
 Cat No :RP10795  
 Lot No :E8357721F  
 Time Processed :20:28:55  
 Month-Day-Year Processed :07/05/2021


Pump A : 0.065% trifluoroacetic in 100% water (v/v)  
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)  
 Total Flow:1 ml/min  
 Wavelength:220 nm

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm



Peak Table

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	15.984	28355	2369	0.170
2	16.523	16453483	2047526	98.472
3	16.913	227008	11102	1.359
Total		16708846	2060998	100.000