


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 04/12/19
Time Processed : 08:34:43 AM
Injection Volume : 0.2
Sample Name : Jagged-1 (188-204)
Cat No. : RP20331
Lot No. : PE3321904
Theoretical MW : 2107.40
Observed MW : 2107.2

Interface : ESI

Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment

Interface Bias : -3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH

Equipment : GK11010007