


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 02/09/17
Time Processed : 17:15:47
Injection Volume : 0.2
Sample Name : IRBP (1-20), human
Cat. No. : RP20269
Lot. No. : PE4821702
Theoretical MW : 2194.60
Observed MW : 2194.40

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH