


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 11/21/17  
Time Processed : 18:41:20  
Injection Volume : 0.2  
Sample Name : DYKDDDDDKPeptide  
Cat No. : RP10586  
Lot No. : PE5341711  
Theoretical MW : 1012.97  
Observed MW : 1013.00

## Interface

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250 °C  
Block Temp : 200 °C

## Interface Bias

Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH