

Sample Name :Glucagon-Like Peptide (GLP) I (7-37)

Cat. No. : RP12738 Lot. No. : PE5721808

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Pump B Conc.	5
25.00	Pumps	Pump B Conc.	65
25.01	Pumps	Pump B Conc.	95
27.00	Pumps	Pump B Conc.	95
27.01	Pumps	Pump B Conc.	5
32.00	Pumps	Pump B Conc.	5
32.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	6.678	27064	140	0.948
2	11.458	1488	177	0.052
3	12.150	10925	1790	0.383
4	15.879	59613	4427	2.087
5	16.322	2736425	266372	95.813
6	16.887	7920	766	0.277
7	17.380	10771	1885	0.377
8	18.325	1806	241	0.063
Total		2856011	275799	100.000