

Sample Name :Melittin,honey bee  
 Cat. No. :RP20415  
 Lot. No. :PE2591707

Pump A : 0.065% trifluoroacetic in 100% water (v/v)  
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)  
 Total Flow:1 ml/min  
 Wavelength:220 nm


Time	Unit	Command	Value	Comment
0.01	Pumps	Pump A B.Conc	5	
25.00	Pumps	Pump A B.Conc	65	
25.01	Pumps	Pump A B.Conc	95	
31.00	Pumps	Pump A B.Conc	95	
31.01	Pumps	Pump A B.Conc	5	
40.00	Pumps	Pump A B.Conc	5	
40.01	Controller	Stop		

<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm

Chromatogram



1 Det.A Ch1 / 220nm

Peak Table

Detector A Ch1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	15.986	4729	532	0.127
2	16.514	26247	1763	0.707
3	16.889	13334	1264	0.359
4	17.122	14005	1045	0.377
5	17.551	3571425	415584	96.250
6	17.878	79728	7863	2.149
7	18.150	1102	409	0.030
Total		3710570	428460	100.000