


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 11/27/16
Time Processed : 15:37:01
Injection Volume : 0.2
Sample Name : Adrenomedullin (AM) (1-52), human
Cat No. : RPI1288
Lot No. : PE1341611
Theoretical MW : 6028.76
Observed MW : 6028.75

Interface

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH