


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 12/11/17
Time Processed : 14:00:03
Injection Volume : 0.2
Sample Name : RGD
Cat. No. : RP20297
Lot. No. : PE5501712
Theoretical MW : 1025.06
Observed MW : 1025.30

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias
Drying Gas Flow
T.Flow
B.conc

: +4.5 kV
: 5 L/min
: 0.2 ml/min
: 50% H2O/50% MeOH