


# Mass Spectrum



## Sample Information

Acquired by : Gary  
Month-Day Processed : 05/23/17  
Time Processed : 10:09:22  
Injection Volume : 0.2  
Sample Name : Calcitonin Gene Related Peptide (CGRP) (8-37), rat  
Cat No. : RP11090  
Lot No. : PE6651705  
Theoretical MW : 3127.60  
Observed MW : 3127.50

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50%H<sub>2</sub>O/50%MeOH