


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 02/06/20  
Time Processed : 06:03:02 PM  
Injection Volume : 0.3  
Sample Name : Calcitonin Gene Related Peptide (CGRP), rat  
Cat No. : RP11095  
Lot No. : E6099420B  
Theoretical MW : 3806.25  
Observed MW : 3806.4

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H<sub>2</sub>O/50%MeOH