


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 01/20/18
Time Processed : 13:14:08
Injection Volume : 0.2
Sample Name : β -Amyloid (1-40)
Cat. No. : RP10004
Lot. No. : PE1651801
Theoretical MW : 4329.82
Observed MW : 4330.60

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias
Drying Gas Flow
T.Flow
B.conc

: +4.5 kV
: 5 L/min
: 0.2 ml/min
: 50%H2O/50%MeOH