


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 10/12/18
Time Processed : 9:27:56
Injection Volume : 0.2
Sample Name : β -Amyloid (25-35)
Cat No. : RP10008
Lot No. : PE2441810
Theoretical MW : 1060.27
Observed MW : 1060.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH