

Sample Name : β -Amyloid (1-42), human

Cat No.: RP20527

Lot No.: PE3721805

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow: 1 ml/min

Wavelength: 220 nm


Time	Unit	Command	Value	Comment
0.01	Pumps	Pump B Conc.	10	
25.00	Pumps	Pump B Conc.	100	
25.01	Pumps	Pump B Conc.	95	
31.00	Pumps	Pump B Conc.	95	
31.01	Pumps	Pump B Conc.	10	
40.00	Pumps	Pump B Conc.	10	
40.01	Controller	Stop		

<<Column Performance>>

<Detector A>

Column : UniPS 10-300 4.6 x 250 mm

Chromatogram



1 Det.A Ch1 / 220nm

Peak Table

Detector A Ch1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	7.510	11964	1049	0.041
2	12.077	194886	21598	0.674
3	12.241	235644	22213	0.815
4	12.664	27617185	1675595	95.469
5	13.644	688305	35595	2.379
6	14.481	131538	14603	0.455
7	14.863	38650	4350	0.134
8	20.537	7526	497	0.026
9	22.231	2290	153	0.008
Total		28927988	1775652	100.000