


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 11/03/17
Time Processed : 15:16:20
Injection Volume : 5
Sample Name : C-Peptide (3-33), human
Cat No. : RP11116
Lot No. : PE0281710
Theoretical MW : 3020.26
Observed MW : 3020.30

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH