


Mass Spectrum



Sample Information
 Month-Day Processed : 12/26/20
 Time Processed : 9:35:38 PM
 Injection Volume : 0.4
 Sample Name : Neuropeptide Y, human
 Cat No : RP10492
 Lot No : E4283620L
 Theoretical MW : 4271.68
 Observed MW : 4270.5

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment : GK11010007
 Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50%H₂O/50%MeOH