


# Mass Spectrum



Sample Information  
Month-Day Processed : 06/23/21  
Time Processed : 2:57:30 PM  
Injection Volume : 0.4  
Sample Name : LL-37, scrambled  
Cat No : RP20332  
Lot No : E8685221F  
Theoretical MW : 4493.30  
Observed MW : 4493.3

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50%H2O/50%MeOH