


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 05/11/18
Time Processed : 10:50:05
Injection Volume : 0.2
Sample Name : CGRP (8-37), rat
Cat. No. : RP11090
Lot. No. : PE5351805
Theoretical MW : 3127.60
Observed MW : 3127.40

Interface

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH