

Sample Name : CEF20
 Cat No: RP11781
 Lot No : E7892720L
 Time Processed : 12:41:07
 Month-Day-Year Processed : 12/23/2020

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow: 1 ml/min
 Wavelength: 220 nm

<<LC Time Program>>


Time	Module	Command	Value
0.01	Pumps	B.Conc	0
25.00	Pumps	B.Conc	60
25.01	Pumps	B.Conc	95
31.00	Pumps	B.Conc	95
31.01	Pumps	B.Conc	0
40.00	Pumps	B.Conc	0
40.00	Controller	Stop	

<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm
 Equipment: ZJ19010140

<Chromatogram>



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	12.773	98956	10516	0.962
2	13.159	9831286	1427273	95.528
3	13.379	60477	12189	0.588
4	13.579	54594	11510	0.530
5	13.642	127551	19239	1.239
6	13.794	21931	4253	0.213
7	13.953	96772	15886	0.940
Total		10291567	1500866	100.000