


# Mass Spectrum



<p><b>Sample Information</b>          Month-Day Processed : 05/28/21          Time Processed : 7:22:22          Injection Volume : 0.4          Sample Name : LL37 (Human)          Cat. No. : RP13323          Lot. No. : E4434921E          Theoretical MW : 4493.28          Observed MW : 4492.8</p>	<p><b>Interface</b> : ESI          Nebulizing Gas Flow : 1.5L/min          CDL Temp : 250          Block Temp : 200</p> <p><b>Equipment</b> : GK11010007          Interface Bias : +4.5 kV          Drying Gas Flow : 5 L/min          T.Flow : 0.2 ml/min          B.conc : 50% H<sub>2</sub>O/50% MeOH</p>
--	--