


Mass Spectrum



Sample Information

Month-Day Processed : 11/21/20

Time Processed : 21:04:33

Injection Volume : 0.3

Sample Name : Luteinizing Hormone Releasing Hormone (LH-RH), human

Sample ID : U7754FK160-1

Theoretical MW : 1182.29

Observed MW : 1181.8

Interface

Interface : ESI

Nebulizing Gas Flow : 1.5L/min

CDL Temp : 250

Block Temp : 200

Equipment

Equipment : GK11010007

Interface Bias : +4.5 kV

Drying Gas Flow : 5 L/min

T.Flow : 0.2 ml/min

B.conc : 50% H₂O/50% MeOH