


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 07/11/20
Time Processed : 22:05:54
Injection Volume : 0.3
Sample Name : OVA G4 peptide
Cat No. : RP20398
Lot No. : E9004720G
Theoretical MW : 906.08
Observed MW : 906.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH