


Mass Spectrum



Sample Information

Acquired by : Gary
 Month-Day Processed : 11/28/17
 Time Processed : 19:25:52
 Injection Volume : 0.2
 Sample Name : C-Peptide (3-33), human
 Cat. No. : RP11116
 Lot. No. : PE5321711
 Theoretical MW : 3020.26
 Observed MW : 3020.40

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250°C
 Block Temp : 200°C

Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50% H2O/50% MeOH