


Mass Spectrum



Sample Information
Month-Day/Processed : 06/18/21
Time Processed : 17:45:34
Injection Volume : 0.4
Sample Name : hgp100(25-33)
Cat No. : RE20344
Lot No. : E4054221F
Theoretical MW : 1155.30
Observed MW : 1155.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH