


Mass Spectrum



Sample Information

Month-Day/Processed : 12/22/20
Time Processed : 10:46:22 PM
Injection Volume : 0.4
Sample Name : hgp100(25-33)
Cat No. : RP20344
Lot No. : E4921420L
Theoretical MW : 1155.30
Observed MW : 1155.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH